


# DURATWIX® JUNIOR NEO / NCF

- ▶ Tracheostomy tube for neonates, infants and toddlers (paediatrics), available without CUFF (REF 17432) and with CUFF (REF 17434)
- ▶ Adaptive, supportive and lightweight material
- ▶ Free of latex and DEHP
- ▶ Soft, flexible neck flange with imprint of the essential specifications of the cannula
- ▶ Rounded tube tip to prevent irritation of the mucous membranes in the trachea
- ▶ With X-ray contrast strip
- ▶ Set includes a broad, stable tube holder, an obturator and cannula passport



## Order information:

**DURATWIX® JUNIOR NEO**

Size 2,5-4



**REF 17432**

**DURATWIX® JUNIOR NCF**

Size 2,5-4



**REF 17434**

When ordering, please state the size after the item number!

Size	O.D. mm	I.D. mm	Angle ca.	Cuff (Ø) mm	Length over outer curve mm	Length over middle curve mm
2,5	4	2,5	120°	8	32	30
3	4,5	3	120°	9	33	31
3,5	5,1	3,5	120°	10	36	33
4	5,8	4	120°	12	38	35

O.D. = Outer diameter · I.D. = Inner diameter