

# InHealth 22.5 Fr Indwelling Voice Prosthesis

sterile

REF VP22xx-LEF



- ▶ The InHealth 22.5 FR Indwelling Voice Prosthesis is an indwelling prosthesis with a low-pressure valve
- ▶ Packed sterile
- ▶ Silver oxide-containing silicone in the valve section reduces candida-related early valve replacement
- ▶ Features a stabilising titanium ring in the shaft
- ▶ With large oesophageal flange, which may be indicated in the presence of periprosthetic leaks
- ▶ The voice prosthesis is placed using the gel capsule technique
- ▶ Diameter 22.5 French
- ▶ Available in sizes 6, 8, 10 and 12

## SCOPE OF DELIVERY

### InHealth 22.5 Fr Indwelling Voice Prosthesis

REF VP22xx-LEF

- 1 x Voice prosthesis (sterile)
- 1 x Flange introducer (sterile)
- 1 x Insertion stick (sterile)
- 2 x Gel capsules
- 1 x Manual

Size	AD ø French	Shaft-thickness mm	Length mm	REF
<b>06</b>	22.5	7,5	6,0	<b>VP2206-LEF</b>
<b>08</b>	22.5	7,5	8,0	<b>VP2208-LEF</b>
<b>10</b>	22.5	7,5	10,0	<b>VP2210-LEF</b>
<b>12</b>	22.5	7,5	12,0	<b>VP2212-LEF</b>

Andreas Fahl Medizintechnik-Vertrieb GmbH is the exclusive distributor for InHealth Technologies® in selected European countries. InHealth Technologies® is a business unit of Freudenberg Medical LLC.