


# LARYVOX<sup>®</sup> DISTANCE HOLDER

- ▶ To be pushed onto the neck flange of a tracheostomy tube:  
ensures the necessary spacing between the cannula opening and, e.g. a stoma neckerchief (e.g. TRACHEOTEX<sup>®</sup> Scarf) or a stoma protection bib (e.g. TRACHEOTEX<sup>®</sup> BIB)
- ▶ Thus an obstruction of the opening of the tracheostomy tube is avoided and unimpaired breathing is possible
- ▶ Made of stainless steel



Single patient use



Consult instructions for use



CE marking

## Order information:

**LARYVOX<sup>®</sup> Distance holder**

**REF 34000**