



# LARYTAPE® / LARYTAPE® TOUCH

- ▶ Adhesive foam filter for laryngectomized patients to provide cover and protection for the tracheostoma
- ▶ It further humidifies, filters and warms the inhaled air to a certain degree
- ▶ Its rounded shape assures an excellent fit on the tracheostoma
- ▶ The circular uninterrupted adhesive stripe assures a tight fit to the tracheostoma without any gaps; thus the foam of LARYTAPE® assumes the function of an HME with mechanical protection from dust particles or insects as well as to retain the exhaled humidity and warmth
- ▶ The adhesive strength may further be enhanced by use of the OPTIGARD® protective dressing wipe (cf separate data sheet)
- ▶ Comes in different sizes (S, L and XL) and in an inconspicuous beige skin colour
- ▶ **Special model: LARYTAPE® TOUCH with push button to facilitate finger occlusion of the tracheostoma when speaking with a voice prosthesis**



## Order information:

LARYTAPE® S, 4 x 5.5 cm, Pack of 10

REF 45960

LARYTAPE® L, 4.8 x 6.7 cm, Pack of 10

REF 45962

LARYTAPE® XL, 6 x 8.5 cm, Pack of 10

REF 45963

LARYTAPE® TOUCH S, 4 x 5.5 cm, Pack of 10

REF 45970

LARYTAPE® TOUCH L, 4.8 x 6.7 cm, Pack of 10

REF 45972

LARYTAPE® TOUCH XL, 6 x 8.5 cm, Pack of 10

REF 45973

- ⊗ Do not reuse
- ☂ Keep dry
- ☀ Keep away from sunlight
- 📖 Consult instructions for use
- CE CE marking