



TRACHLINE® UNI CUFF-SUCTION

1 IC + 2 ICU

- ▶ Tracheostomy tube with CUFF and integrated suction tube for long-term use, max. service life up to 29 days
- ▶ Soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ Rounded tube tip to prevent irritation of the mucous membranes
- ▶ Bayonet locks enable simple, secure attachment of the inner cannula and accessories
- ▶ Set consists of an inner cannula (IC), two inner cannulas with 15 mm connector (ICU), a 15 mm connector (UNI), a cough cap, a locking ring, an adjustable tube holder, an obturator as well as a cannula passport



Order information:

TRACHLINE® UNI CUFF-SUCTION
1 IC + 2 ICU

Size 7-9



REF 13802

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	I.D. O.C. mm	I.D. I.C. mm	Angel ca.	Cuff (Ø) mm	Length mm
7	13,1	8.5	7.0	105°	22	72
8	14,5	9.5	8.0	100°	26	80
9	15,4	10.5	9.0	95°	30	82

O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · O.C. = Outer cannula · I.C. = Inner cannula · Length = Length over middle curve