

LARYVOX® TWEEZERS

Encrustation tweezers

- ▶ To remove encrustations from the tracheostoma:
The encrustations of tough secretions can be quickly and safely gripped and removed from the tracheostoma
- ▶ The rounded tips enable gentle operation and minimize the risk of injury
- ▶ The encrustation tweezers are angled so that an unimpaired view of the tracheostoma in front of a mirror is possible
- ▶ Made of stainless steel



Consult instructions for use



CE marking

Order information:

LARYVOX® Tweezers

REF 35000