

LARYVOX® EMERGENCY-SET TRACHEOSTOMY

Complete set for intensive care

- ▶ 1x LARYVOX® SOS safe-mask
Emergency ventilation mask for laryngectomees
- ▶ 1x LARYVOX® Tracheospreizer (tracheal dilator)
Using the LARYVOX® Tracheospreizer, the cannula can be inserted safely and as quickly as possible
- ▶ 1x HUMIDOTRACH®
HME filter with oxygen supply port and 15 mm standard connector; fits to all tracheostomy tubes with 15 mm connector (UNI)
- ▶ 1x OPTIFLAUSCH® K
“K” variant: with adherent fastener assures the secure fixation of tracheostomy tubes
- ▶ 2x 2-KAM® SLIT, sterile Kompresse
Tracheal compresses consist of two layers pad the cannula softly and saely against the tracheostoma and the skin and provides a firm sit
- ▶ 3x OPTIFLUID® Stoma oil towel
Stoma oil towels with colourless neutral oil to lubricate the tracheostomy tube before insertion
- ▶ 3x OPTIFAHL® Stoma cleaning wipes
Cleaning wipe for the cleaning and care of the skin around the tracheostoma
- ▶ 1x Oxygen connection tube
- ▶ 1x Single-use syringe / 10ml



For all components of this kit separate and detailed data sheets are available.

Order information:

LARYVOX® EMERGENCY-SET Tracheostomy

REF 48008