


















PRODUCT RANGE LARYVOX® TAPE

REF	Product Name	Picture	Material Properties
48100	LARYVOX® TAPE STANDARD ROUND		<ul style="list-style-type: none"> • Good adhesive strength • Transparent appearance • Water-repellent effect • Suitable in situations with unproblematic tracheostoma anatomy
48200	LARYVOX® TAPE STANDARD OVAL		
48300	LARYVOX® TAPE STANDARD XL OVAL		
48120	LARYVOX® TAPE FLEXIBLE ROUND		<ul style="list-style-type: none"> • High adhesive strength • Suitable for problematic tracheostoma situations • Suitable for normal physical activity
48220	LARYVOX® TAPE FLEXIBLE OVAL		
48320	LARYVOX® TAPE FLEXIBLE XL OVAL		
48130	LARYVOX® TAPE HYDROSOFT ROUND		<ul style="list-style-type: none"> • Particularly skin-friendly • High fixation comfort • Ability to enclose fluids/secretions
48230	LARYVOX® TAPE HYDROSOFT OVAL		
48330	LARYVOX® TAPE HYDROSOFT XL OVAL		

REF	Product Name	Picture	Material Properties
48140	LARYVOX® TAPE COMFORT ROUND		<ul style="list-style-type: none"> • Highest possible adhesive properties • Adapts comfortably to movements • High flexibility and long wearing period • Optimal for physical activity
48240	LARYVOX® TAPE COMFORT OVAL		
48340	LARYVOX® TAPE COMFORT XL OVAL		
48050	LARYVOX® TAPE HYPOALLERGEN ROUND		<ul style="list-style-type: none"> • Particularly for sensitive skin • Reduces allergic skin reactions • Adapts comfortably to movements • Very good adhesive properties • High flexibility and long wearing period
48150	LARYVOX® TAPE HYPOALLERGEN OVAL		
48250	LARYVOX® TAPE HYPOALLERGEN XL OVAL		
48400	LARYVOX® TAPE EXTRA FINE		<ul style="list-style-type: none"> • Shape: rectangular • Very good adhesive properties • Transparent and extremely thin base plate • Water-repellent material • Covers a large skin area around the tracheostoma • Adapts optimally to the movements of the neck region • Discreet appearance
48500	LARYVOX® TAPE CONVEX		<ul style="list-style-type: none"> • Suitable for precise positioning in a deep or deepened tracheostoma • The base plate as honeycomb core offers a strong overall structure • Reduces the occurrence of allergic skin reactions • Transparent appearance • Water-repellent material • Very good adhesive properties

ANDREAS FAHL

MEDIZINTECHNIK-VERTRIEB GMBH

August-Horch-Straße 4a · 51149 Köln · Phone +49(0)22 03/29 80-0
 Fax +49(0)2203/2980-100 · Germany · mail vertrieb@fahl.de · www.fahl.de