


OPTIFLAUSCH® K JUNIOR PED

REF 32556



REF 32558



- ▶ Assures the secure fixation of tracheostomy tubes
 - “K PED” variant: with adherent fastener, for pediatric use
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Broad and soft for high wearing comfort
- ▶ Approx. 4 mm thick and 2.5 cm or 1.7 cm wide.
- ▶ The tube holder is wrapped around the neck and the adherent fastener enables it to be easily and quickly adjusted to the desired neck size (2 parts, adjustable).
- ▶ Colour: animal design
- ▶ Individually packed

SCOPE OF DELIVERY

- ▶ Manual

OPTIFLAUSCH® K JUNIOR PED **REF 32556**
1 x Tracheostomy Tube Holder

or

OPTIFLAUSCH® SLIM K JUNIOR PED **REF 32558**
1 x Tracheostomy Tube Holder

Product	Lenght	Width	REF
OPTIFLAUSCH® K JUNIOR PED	16 - 25 cm	2,5 cm	32556
OPTIFLAUSCH® SLIM K JUNIOR PED	16 - 30 cm	1,7 cm	32558

