

DURATWIX® UNI-VARIO

1 ICU + 1 ICV

- ▶ Tracheostomy tube with two inner cannulas (IC), both with 15 mm connector, one with rotating 15 mm connector (ICV) to buffer torque and tensions
- ▶ Soft double-swivel ROTATWIX® neckflange, rotating in 360° for enhanced neck mobility
- ▶ Medical-grade polyurethane, soft, atraumatic, thermosensible
- ▶ Modular: 6 different types of inner cannulas
- ▶ Set includes broad and adjustable tube holder (OPTIFLAUSCH®), obturator and product ID card
- ▶ Size 7-10 available in two lengths
- ▶ Colour coding
 Inner cannula: fenestrated/unfenestrated, standard length/short
 Imprint outer cannula neckflange: sieve-fenestrated/non-fenestrated



STERILE EO Sterilized using ethylene oxid

Do not resterilize

Keep away from sunlight

Single patient use

Consult instructions for use

Rx only CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

CE 0482 CE marking identification number of the notified body


Product made in Germany

Order information:

DURATWIX® UNI-VARIO 1 ICU + 1 ICV, SIZE 6-10

DURATWIX® UNI-VARIO short 1 ICU + 1 ICV, SIZE 7-10

When ordering, please state the size after the item number!



REF 19211

REF 19212

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm standard	REF standard	Length mm short	REF short
6	11.9	10.3	6.0	82.0	19211-06	-	-
7	11.9	10.8	7.0	82.0	19211-07	72.0	19212-07
8	12.5	11.4	8.0	85.0	19211-08	75.0	19212-08
9	13.7	12.5	9.0	88.0	19211-09	78.0	19212-09
10	15.0	13.8	10.0	91.0	19211-10	81.0	19212-10

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve