

TRACHFLOW® METRIC S + L

sterile

REF 68630/68730



REF 60700




- ▶ Sterile, atraumatic tracheostoma catheters for particularly gentle suctioning of tracheal secretions
- ▶ With rounded tip to minimize the risk of tracheal injury
- ▶ Made of soft, flexible and transparent plastic
- ▶ With central opening and four lateral openings ("eyes") at the tip to avoid adherence to the trachea
- ▶ With grading in 1 cm increments to facilitate control of the insertion depth
- ▶ Wide range of sizes (CH 10-16)
- ▶ Colour coded connectors Available in two different lengths: 25 cm and 50 cm (the short variant is particularly suited for tracheal suctioning)
- ▶ To be used with the TRACHFLOW® Fingertip, to control the suction power. The Fingertip is connected between the suction device and the sterile stoma catheter. One size of the TRACHFLOW® Fingertip fits all sizes of stoma catheters

SCOPE OF DELIVERY

- ▶ Manual
- | | |
|---|------------------|
| TRACHFLOW® METRIC S | REF 68630 |
| 1 x Suction catheter, sterile, atraumatic | |
| or | |
| TRACHFLOW® METRIC L | REF 68730 |
| 1 x Suction catheter, sterile, atraumatic | |
| Optional accessories: | |
| TRACHFLOW® FINGERTIP | REF 60700 |

INFORMATION



Article	Size	Length	Shore	Color	Unit	REF
TRACHFLOW® METRIC S	CH 10	25 cm	A64	Black	30	68630-10
TRACHFLOW® METRIC S	CH 12	25 cm	A64	White	30	68630-12
TRACHFLOW® METRIC S	CH 14	25 cm	A64	Green	30	68630-14
TRACHFLOW® METRIC S	CH 16	25 cm	A64	Orange	30	68630-16
TRACHFLOW® METRIC L	CH 10	50 cm	A64	Black	30	68730-10
TRACHFLOW® METRIC L	CH 12	50 cm	A64	White	30	68730-12
TRACHFLOW® METRIC L	CH 14	50 cm	A64	Green	30	68730-14
TRACHFLOW® METRIC L	CH 16	50 cm	A64	Orange	30	68730-16

CE 0482 STERILE EO