

OPTISWAB® SET PURE / OPTISWAB® SET PURE XL

REF 31025



REF 31030



- ▶ Both sets for cleaning the tracheostomy tubes include:
 - **OPTIBRUSH® CONT XL** Cleaning tub (800 ml) (REF 31300) with large diameter to enable cleaning of the tracheostomy tubes. The cleaning tub contains a removable sieve insert which allows moving the tubes in the solution.
 - **OPTICIT®** liquid disinfectant (concentrate /250 ml) (REF 31180): based on glutardialdehyde, a proven disinfectant active ingredient for the disinfection of medical products. Extremely effective against endospores, many microorganisms and viruses, incl. MRSA, HBV/HIV/HCV). Easy to use by diluting with water to create the ready-to-use solution.
- ▶ **OPTISWAB® SET PURE (REF 31025)** contains:
 - **OPTIBRUSH® SWAB** cleaning swabs (REF 31910) for the gentle cleaning of tracheostomy tubes. Flexible plastic swab with absorbent polyurethane foam cap. No fluff or fibres. This reduces the risk of tracheal irritation. Particularly suitable for tracheostomy tubes made of soft, sensitive plastic, especially for paediatric tracheostomy tubes.
- ▶ **OPTISWAB® SET PURE XL (REF 31020)** contains:
 - **OPTIBRUSH® SWAB XL (REF 31920)** cleaning swabs for the gentle cleaning of tracheostomy tubes. Flexible plastic swab with absorbent polyurethane foam cap. No fluff or fibres. This reduces the risk of tracheal irritation. Particularly suitable for tracheostomy tubes made of soft, sensitive plastic.

SCOPE OF DELIVERY

- ▶ Instruction of use

OPTISWAB® SET PURE REF 31025

1 x OPTIBRUSH® CONT XL cleaning tub
 1 x OPTICIT® disinfectant concentrate
 4 x OPTIBRUSH® SWAB cleaning swabs

or

OPTISWAB® SET PURE XL REF 31030

1 x OPTIBRUSH® CONT XL cleaning tub
 1 x OPTICIT® disinfectant concentrate
 4 x OPTIBRUSH® SWAB cleaning swabs

All set components are also available individually.

Product

OPTISWAB® SET PURE

REF

31025

OPTISWAB® SET PURE XL

31030