

**SENSOFOAM®/ SENSOFOAM® PAD**

sterile


REF 30860



REF 30861



REF 30862



- ▶ High-quality foam compress with high absorption capacity
- ▶ Made of a 3-layer foam material to prevent adhesion to the skin and to absorb secretions
- ▶ Due to its reduced adhesion to the skin the compress may be removed easily
- ▶ Transmissible to air and humidity to prevent skin maceration and inflammation
- ▶ With slit to allow easy exchange of the compress while the tracheostomy tube remains in situ
- ▶ Stabilizes tight seal of the tracheostoma and a tight fit of the tracheostomy tube in the trachea
- ▶ Sterile to allow immediate post-operative use

**SCOPE OF DELIVERY**

- ▶ Manual

**SENSOFOAM®** **REF 30860**  
10 x tracheal compresses, sterile

or  
**SENSOFOAM®** **REF 30861**  
10 x tracheal compresses, sterile, small

or  
**SENSOFOAM® PAD** **REF 30862**  
10 x tracheal compresses, sterile

Product	Size	Unit	REF
<b>SENSOFOAM®</b>	10,0 x 9,0 x 0,5 cm	10	<b>30860</b>
<b>SENSOFOAM®</b>	6,5 x 6,5 x 0,5 cm	10	<b>30861</b>
<b>SENSOFOAM® PAD</b>	7,8 x 3,8 cm	10	<b>30862</b>

CE 0482 STERILE