



DURAPUR® LINGO-PHON

1 PHON ICF



- ▶ Tracheostomy tube for long-term use, can be used up to 6 months
- ▶ Outer cannula with sieve-fenestration, with inner cannula fenestrated with silver speaking valve (PHON ICF) and flat profile if no connection to any peripheral device (ventilation equipment, HME or speaking valve) is desired
- ▶ The DURAPUR® tracheostomy tubes are made of thermosensitive medical synthetics that develop their optimal product properties at body temperature
- ▶ Rounded tube tip to prevent irritation of the mucous membranes in the trachea
- ▶ Transparent, inconspicuous design
- ▶ Tracheostomy tube available in size 7-12
- ▶ Including broad adjustable tube holder (NECKFIX®) and a product ID card
- ▶ Inner cannula (PHON-ICF) also available separately (REF 23808), size to be selected according to the size of the DURAPUR outer cannula

Order information:

DURAPUR® LINGO-PHON 1 ICF SIZE 7-12



REF 23032

DURAPUR® PHON ICF SIZE 7-12



REF 23808

When ordering, please state the size after the item number!

- Keep away from sunlight
- Single patient use
- Consult instructions for use
- CE marking identification number of the notified body

Size	O.D. neck flange mm	O.D. tip mm	I.D. O.C. mm	I.D. I.C. mm	Length mm standard	Length mm extra-lang (XL)
7	10.6	8.7	7.2	5.6	65	-
8	12.5	10.4	8.2	6.1	70	105
9	13.6	11.3	9.3	7.1	75	110
10	14.6	12.4	10.0	8.0	85	115
11	15.3	13.3	11.1	9.0	87	-
12	16.5	14.3	12.3	10.2	90	-

O.D. neck flange = Outer diameter behind the neck flange · O.D. tip = Outer diameter of the outer cannula at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula · O.C. = Outer cannula · I.C. = Inner cannula
Length = Length over the outer curve