

# HUMIDOPHONE® PLUS / HUMIDOPHONE® PLUS O<sub>2</sub>

REF 46487 CE



REF 46489 CE 0482



- ▶ Speaking Valve with integrated filter
- ▶ Membrane closed in basic position ("closed-position valve") – opens during inspiration
- ▶ Advantage: No enhanced exhalation pressure required to close the valve and to facilitate speaking. Less fatigue for the patient, better speaking quality
- ▶ Low profile, lightweight
- ▶ Open-pore polyurethane filter
- ▶ Filter in proximal position (towards the patient)
- ▶ With multi-adapter to fit both 22 mm combi-adapter devices (KOMBI) and 15 mm connector devices (UNI)
- ▶ **HUMIDOPHONE® PLUS O<sub>2</sub>** (REF 46489) with oxygen supply port (plug-in connection 5 mm)

## SCOPE OF DELIVERY

- ▶ Instruction for use

**HUMIDOPHONE® PLUS**                      **REF 46487**  
30 x Speaking Valve

or  
**HUMIDOPHONE® PLUS O<sub>2</sub>**                      **REF 46489**  
30 x Speaking Valve

The HUMIDOPHONE® speaking valve must be only used by patients with a partially or fully preserved larynx!

Product	Unit	REF
<b>HUMIDOPHONE® PLUS</b>	30 pieces, packed by 5 pcs.	<b>46487</b>
<b>HUMIDOPHONE® PLUS O<sub>2</sub></b>	30 pieces, packed by 5 pcs.	<b>46489</b>