


SENSOTRACH®


REF 30520




REF 30609



REF 30662



REF 30780



- ▶ High-quality tracheal compresses with anatomical rounded shape for maximum wearing comfort
- ▶ **UNO** variants single-layered, for patients with little production of secretion
- ▶ **DUO** variants with two layers, for patients with moderate to average production of secretion: the layer close to the skin – marked by a green thread - effectively absorbs secretions and transfers them to the peripheral layer. This keeps the skin dry and warm
- ▶ **3-PLUS** variants with three layers, for patients with normal to strong production of secretion: the layer towards the skin effectively absorbs secretions and transfers them to the central and peripheral layer. This keeps the skin dry and warm. The peripheral layer has a moisture-repellent to protect the clothing
- ▶ **DUO ALU** variants are a high quality, especially absorbent and dimensionally stable tracheal compresses, which have an aluminium coating on one side. The metallic coating prevents sticking to the skin, due to blood and secretions.
- ▶ Standard variants without slit
- ▶ **SLIT** variants with slit for exchange while the tracheostomy tube is in situ

SCOPE OF DELIVERY

- ▶ Manual

SENSOTRACH® UNO REF 305xx
10 x tracheal compresses

or
SENSOTRACH® DUO REF 3060x
10 x tracheal compresses

or
SENSOTRACH® DUO ALU REF 3066x
10 x tracheal compresses

or
SENSOTRACH® 3-PLUS REF 307xx
10 x tracheal compresses

Product	Size	Unit	REF
SENSOTRACH® UNO	9,0 x 9,8 cm	10	30520
SENSOTRACH® UNO SLIT	9,0 x 9,8 cm	10	30521
SENSOTRACH® DUO	9,0 x 9,8 cm	10	30608
SENSOTRACH® DUO SLIT	9,0 x 9,8 cm	10	30609
SENSOTRACH® DUO ALU	9,0 x 9,8 cm	10	30662
SENSOTRACH® DUO ALU SLIT	9,0 x 9,8 cm	10	30663
SENSOTRACH® 3-PLUS	9,0 x 9,8 cm	10	30770
SENSOTRACH® 3-PLUS SLIT	9,0 x 9,8 cm	10	30780

