

# SUCTION CATHETER IDEAL TIP

REF 019055x1 / REF 019068x1



REF 60700 (optional)



- ▶ Sterile tracheostoma catheters for suctioning tracheal secretions
- ▶ Made of soft, flexible and transparent plastic with rounded tip
- ▶ With central opening and two lateral openings ("eyes") at the tip to avoid adherence to the trachea
- ▶ Shore hardness A74
- ▶ Available in two different lengths: 25 cm and 50 cm
- ▶ Wide range of sizes (CH 06-18)
- ▶ Colour coded connectors
- ▶ To be used with the TRACHFLOW® Fingertip, to control the suction power. The Fingertip is connected between the suction device and the sterile stoma catheter. One size of the TRACHFLOW® Fingertip fits all sizes of suction catheters
- ▶ DEHP and latex free

## SCOPE OF DELIVERY

**SUCTION CATHETER IDEAL TIP, 50 CM** REF 019055X1

150 x Suction catheter, sterile

or

**SUCTION CATHETER IDEAL TIP, 25 CM** REF 019068X1


100 x Suction catheter, sterile

When ordering, please state the size after the item number!

Optional available:

**TRACHFLOW® FINGERTIP** REF 60700  
(unsterile)

## INFORMATION



Article	Size	Length	Shore	Unit	REF
SUCTION CATHETER IDEAL TIP, 50 CM	CH 06	50 cm	A74	150	01905511
SUCTION CATHETER IDEAL TIP, 50 CM	CH 08	50 cm	A74	150	01905521
SUCTION CATHETER IDEAL TIP, 50 CM	CH 10	50 cm	A74	150	01905531
SUCTION CATHETER IDEAL TIP, 50 CM	CH 12	50 cm	A74	150	01905541
SUCTION CATHETER IDEAL TIP, 50 CM	CH 14	50 cm	A74	150	01905551
SUCTION CATHETER IDEAL TIP, 50 CM	CH 16	50 cm	A74	150	01905561
SUCTION CATHETER IDEAL TIP, 50 CM	CH 18	50 cm	A74	150	01905571
SUCTION CATHETER IDEAL TIP, 25 CM	CH 06	25 cm	A74	100	01906811
SUCTION CATHETER IDEAL TIP, 25 CM	CH 08	25 cm	A74	100	01906821
SUCTION CATHETER IDEAL TIP, 25 CM	CH 10	25 cm	A74	100	01906831
SUCTION CATHETER IDEAL TIP, 25 CM	CH 12	25 cm	A74	100	01906841
SUCTION CATHETER IDEAL TIP, 25 CM	CH 14	25 cm	A74	100	01906851
SUCTION CATHETER IDEAL TIP, 25 CM	CH 16	25 cm	A74	100	01906861
SUCTION CATHETER IDEAL TIP, 25 CM	CH 18	25 cm	A74	100	01906871

CE 0044 STERILE EO (not valid for REF 60700)