

OPTIBRUSH® SET / OPTIBRUSH® SET PLUS

REF 31011



REF 31012



- ▶ Both sets for cleaning the tracheostomy tubes include:
 - OPTIBRUSH® CONT Cleaning Box (REF 31200) with a large diameter to allow cleaning of standard tracheostomy tubes. The Cleaning Box contains a removable sieve insert which allows the tracheostomy tube to be gently agitated in the cleaning solution
 - OPTIBRUSH® CLEAN Cleaning Powder (REF 31110) including measuring spoon for easy and accurate preparation of the cleaning solution; simply dissolve in warm tap water; for a thorough but gentle removal of blood and secretion from tracheostomy tubes
- ▶ **OPTIBRUSH® SET (REF 31011)** contains:
 - OPTIBRUSH® Cleaning Brushes are flexible and have a rounded tip to avoid scratching of the soft tracheostomy tube material
- ▶ **OPTIBRUSH® SET PLUS (REF 31012)** contains:
 - OPTIBRUSH® PLUS (REF 31855) cleaning brushes are flexible and surrounded by a soft, active wiping and absorbent fibre top.
- ▶ The brush size corresponds to the set size.

SCOPE OF DELIVERY

- ▶ Manual

OPTIBRUSH® SET **REF 31011**
 1 x OPTIBRUSH® CONT cleaning tin
 1 x OPTIBRUSH® CLEAN cleaning powder (100 g)
 4 x OPTIBRUSH® cleaning brushes

or

OPTIBRUSH® SET PLUS **REF 31012**
 1 x OPTIBRUSH® CONT cleaning tin
 1 x OPTIBRUSH® CLEAN cleaning powder (100 g)
 4 x OPTIBRUSH® PLUS cleaning brushes

When ordering, please state the size after the item number!

All components of the set are also available individually.

Product	Size	REF
OPTIBRUSH® SET	08	31011-08
OPTIBRUSH® SET	10	31011-10
OPTIBRUSH® SET	12	31011-12
OPTIBRUSH® SET PLUS	08	31012-08
OPTIBRUSH® SET PLUS	10	31012-10
OPTIBRUSH® SET PLUS	12	31012-12