


**SENSOTRACH® DUO SKIN**

sterile

REF 30668



REF 30669



- ▶ High-quality tracheal compresses - beige
- ▶ With anatomical rounded shape for maximum wearing comfort
- ▶ The combination of the anatomical shape and naturally appearing colour ensures that the tracheal compress can barely be seen or felt
- ▶ With two layers: inner layer effectively absorbs secretions and transfers them to the peripheral layer.
- ▶ Keeps the skin dry and warm
- ▶ Available with **SLIT**, for exchange while the tracheostomy tube is in situ: **SENSOTRACH® DUO SKIN SLIT**

**SCOPE OF DELIVERY**

- ▶ Manual

**SENSOTRACH® DUO SKIN** REF 30668  
10 x tracheal compresses, sterile

or  
**SENSOTRACH® DUO SKIN SLIT** REF 30669  
10 x tracheal compresses, sterile

Product	Size	Unit	REF
<b>SENSOTRACH® DUO SKIN</b>	9,0 x 9,8 cm	10	<b>30668</b>
<b>SENSOTRACH® DUO SKIN SLIT</b>	9,0 x 9,8 cm	10	<b>30669</b>

CE 0482 STERILE