

LARYNGOTEC® KOMBI CLIP

- ▶ Tracheostomy tube for long-term use, can be used up to 3-6 months
- ▶ Medical grade silicone: soft, pliable, lightweight
- ▶ Bulging ring/Clip behind the 22 mm adapter allows connection to adhesive base plates
- ▶ Transparent, inconspicuous design
- ▶ With 22 mm adapter (KOMBI) to connect HMEs
- ▶ Fenestrated variant (LARYNGOTEC® KOMBI CLIP LINGO, REF 14820) allows further the connection of a speaking valve
- ▶ Size 7 to 13
- ▶ Each size in three different lengths
- ▶ With (LARYNGOTEC® KOMBI CLIP LINGO, REF 14820) or without (LARYNGOTEC® KOMBI CLIP, REF 14830) sieve-fenestration
- ▶ Custom-made sieve-fenestration possible
- ▶ Including product ID card



Order information:

LARYNGOTEC® KOMBI CLIP, SIZE 7-13



REF 14830

LARYNGOTEC® KOMBI CLIP LINGO, SIZE 7-13



REF 14820

When ordering, please state the size after the item number!

Size	O.D. tip mm	O.D. neck flange mm	I.D. mm	Length mm	Length middle radius mm	REF	REF LINGO
7	10.2	10.2	7.2	44	38	14830-0744	14820-0744
7	10.2	10.2	7.2	62	54	14830-0762	14820-0762
7	10.2	10.2	7.2	65	56	14830-0765	14820-0765
8	11.2	11.2	8.2	44	38	14830-0844	14820-0844
8	11.2	11.2	8.2	62	54	14830-0862	14820-0862
8	11.2	11.2	8.2	70	61	14830-0870	14820-0870
9	12.4	12.4	9.4	44	36	14830-0944	14820-0944
9	12.4	12.4	9.4	62	55	14830-0962	14820-0962
9	12.4	12.4	9.4	75	65	14830-0975	14820-0975
10	13.4	13.4	10.4	44	36	14830-1044	14820-1044
10	13.4	13.4	10.4	62	55	14830-1062	14820-1062
10	13.4	13.4	10.4	85	75	14830-1085	14820-1085
11	14.4	14.4	11.4	44	36	14830-1144	14820-1144
11	14.4	14.4	11.4	62	55	14830-1162	14820-1162
11	14.4	14.4	11.4	87	76	14830-1187	14820-1187
12	15.4	15.4	12.4	44	36	14830-1244	14820-1244
12	15.4	15.4	12.4	62	55	14830-1262	14820-1262
12	15.4	15.4	12.4	90	77	14830-1290	14820-1290
13	16.6	16.6	13.0	44	37	14830-1344	14820-1344
13	16.6	16.6	13.0	62	55	14830-1362	14820-1362
13	16.6	16.6	13.0	90	77	14830-1390	14820-1390

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula
 Length = Length over the outer curve · Length of mean radius = Length on midline, from flange to tip



Sterilized using ethylene oxide



Keep away from sunlight



Single patient use



Consult instructions for use



CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).



CE marking identification number of the notified body

Product made in Germany