


DURAPUR® UNI

1 ICU



- ▶ Tracheostomy tube for long-term use, can be used up to 6 months
- ▶ Inner cannula with 15 mm-UNI-connector (ICU)
- ▶ The DURAPUR® tracheostomy tubes are made of thermosensitive medical synthetics that develop their optimal product properties at body temperature
- ▶ Rounded tube tip to prevent irritation of the mucous membranes in the trachea
- ▶ Transparent, inconspicuous design
- ▶ Available in two different lengths: Standard and Extra-long (XL)
- ▶ Standard length from size 7 to 12, Extra-long (XL) variants from size 8 to 10
- ▶ Including broad adjustable tube holder (NECKFIX®) and a product ID card
- ▶ Inner cannula with standard length (ICU) also available separately (REF 23811), size to be selected according to the size of the DURAPUR outer cannula

Order information:

DURAPUR® UNI 1 ICU	SIZE 7-12		REF 23011
DURAPUR® UNI XL 1 ICU	SIZE 8-10		REF 23211
DURAPUR® ICU	SIZE 7-12		REF 23811

When ordering, please state the size after the item number!

Keep away from sunlight

Single patient use

Consult instructions for use

CE marking identification number of the notified body

Size	O.D. neck flange mm	O.D. tip mm	I.D. O.C. mm	I.D. I.C. mm	Length mm standard	Length mm extra-lang (XL)
7	10.6	8.7	7.2	5.6	65	-
8	12.5	10.4	8.2	6.1	70	105
9	13.6	11.3	9.3	7.1	75	110
10	14.6	12.4	10.0	8.0	85	115
11	15.3	13.3	11.1	9.0	87	-
12	16.5	14.3	12.3	10.2	90	-

O.D. neck flange = Outer diameter behind the neck flange · O.D. tip = Outer diameter of the outer cannula at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula · O.C. = Outer cannula · I.C. = Inner cannula
Length = Length over the outer curve