


TRACHEOSILC® VARIO CUFF



- ▶ Tracheostomy tube with CUFF
- ▶ Also available without CUFF (REF 18820) and with CUFF and integrated suction line (SUCTION) to remove secretions above the CUFF (REF 18840 see separate data sheets)
- ▶ Soft, flexible neck-flange with key tube identification data clearly marked
- ▶ Clear and transparent tube shaft for easy handling during cleaning
- ▶ With rotating 15 mm connector (VARIO) to buffer torque and tension caused by tractive forces from ventilation equipment
- ▶ Rounded cannula tip to prevent irritation of the mucosa in the trachea
- ▶ Set includes tube holder, obturator and product ID card

Order information:

TRACHEOSILC® VARIO CUFF

SIZE 5-9,5



REF 18830

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	O.D. Tip mm	Length mm	REF
5.0	7.3	7.3	57	18830-05
6.0	8.7	8.7	63	18830-06
7.0	10.0	10.0	71	18830-07
7.5	10.7	10.7	73	18830-075
8.0	11.0	11.0	75	18830-08
8.5	11.7	11.7	78	18830-085
9.0	12.3	12.3	80	18830-09
9.5	13.3	13.3	83	18830-095

O.D. neck flange = outer diameter behind the neck flange · O.D. Tip = outer diameter at the tube tip · Length = length over the outer curve