

HUMIDOTWIN®

REF 46460





- ▶ HME filter / Artificial nose with overpressure valve
- ▶ Particularly for nighttime use
- ▶ Multiadapter fits on all 15 mm (UNI) connectors as well as on 22 mm (KOMBI) adapters:
- ▶ Open pore polyurethane foam filter
- ▶ Low profile, light-weight, inconspicuous
- ▶ With overpressure valve to allow unimpaired coughing

SCOPE OF DELIVERY

- ▶ Instruction for use

HUMIDOTWIN®
1 x HME filter cassettes

REF 46460

Product	Unit	REF
HUMIDOTWIN®	30 pieces, individually packed	46460

CE