








# BUCHANAN® BIB

- ▶ To cover and protect the tracheostoma of laryngectomized and tracheostomized patients
- ▶ Particularly suitable for patients with a sensitive skin who do not tolerate adhesive filters
- ▶ A special foam insert ensures effective filtering, warming and humidification of the inspired air
- ▶ Large surface
- ▶ Extra-long neck straps – suits every neck size
- ▶ Two sizes: Large (L) and small (S)
- ▶ Each bib is individually and hygienically packed inside its packaging unit
- ▶ Washable and reusable



	Keep away from sunlight
	Keep dry
	Single patient use
	Consult instructions for use
	CE marking
	Hand wash
	Do not dry in the drum dryer

## Order information:

**BUCHANAN® BIB L**, max. 95 cm, 21 x 19 cm, Pack of 10

**REF 41700**

**BUCHANAN® BIB S**, max. 65 cm, 15 x 12 cm, Pack of 10

**REF 41800**