

DURATWIX® UNI

- ▶ Tracheostomy tube with wide inner lumen to provide maximum volume of air at optimized flow
- ▶ Cannula with 15 mm connector (UNI)
- ▶ Soft double-swivel ROTATWIX® neckflange, rotating in 360° for enhanced neck mobility
- ▶ Medical-grade polyurethane, soft, atraumatic, thermosensible
- ▶ Set includes broad and adjustable OPTIFLAUSCH® tube holder, obturator and product ID card
- ▶ Size 7-10 available in two lengths



- Sterilized using ethylene oxide
- Do not re-sterilize
- Keep away from sunlight
- Single patient use
- Consult instructions for use

Rx only
CAUTION:
 United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

CE 0482 CE marking identification number of the notified body

Product made in Germany

Order information:

DURATWIX® UNI, SIZE 7-10



REF 19201

DURATWIX® UNI short, SIZE 7-10



REF 19202

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm standard	REF standard	Length mm short	REF short
7	11.9	10.8	7.0	82.0	19201-07	72.0	19202-07
8	12.5	11.4	8.0	85.0	19201-08	75.0	19202-08
9	13.7	12.5	9.0	88.0	19201-09	78.0	19202-09
10	15.0	13.8	10.0	91.0	19201-10	81.0	19202-10

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve