



# DURAVENT®

- ▶ Tracheostomy tube for long-term use, can be used up to 3-6 months
- ▶ Soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ Two lengths: Standard length plus extra-long length (XL)
- ▶ Standard length from size 03 to 13, including pediatric sizes with half-mm size 3,5
- ▶ Extra-long (XL) variants from size 05-13
- ▶ Without inner cannula (IC)
- ▶ Custom-made variants possible (length, sieve-fenestration, bending)
- ▶ Including broad adjustable OPTIFLAUSCH® tube holder and a product ID card



## Order information:

**DURAVENT®**, SIZE 3-13

**REF 11000**

**DURAVENT® XL**, SIZE 5-13

**REF 11200**

When ordering, please state the size after the item number!



Size	O.D. neck flange mm	O.D. tip mm	I.D. mm O-cannula	Length mm standard	Length mm extra long (XL)
3	5.2	3.8	2.7	55	-
3.5	6.2	4.8	3.6	55	-
4	7.2	5.8	4.4	55	-
5	8.4	6.6	5.4	55	95
6	9.4	7.8	6.2	60	95
7	10.2	8.8	7.0	65	100
8	11.2	9.8	8.0	70	105
9	12.4	11.2	9.0	75	105
10	13.4	12.2	10.0	85	110
11	14.4	13.2	11.0	87	115
12	15.4	14.2	12.0	90	120
13	16.6	15.2	13.0	90	125

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve

Keep away from sunlight

Single patient use

Consult instructions for use

**CAUTION:**  
United States Federal law  
restricts this device to sale,  
distribution and use by or  
on order of a physician or a  
licensed practitioner (in the  
US only).

Rx only

CE marking identification  
number of the notified body

**Product made in Germany**