


KACLIP®

REF 32000



REF 32005



REF 32011



REF 32015



- ▶ Family of elastic tube holders which assure the secure fixation of tracheostomy tubes
 - “H” variants: with plastic hook
 - “K” variants: with adherent fastener
 - “HM” variants: with metal hook
- ▶ Can be adjusted to neck widths of approx. 30 cm to 50 cm due to its high elasticity
- ▶ Various fixation options: with metal hook, with plastic hook, with adherent fastener
- ▶ Various colour options: white, beige or transparent
- ▶ Individually packed

SCOPE OF DELIVERY

- ▶ Manual
 - KACLIP® HM** REF 3200x
1 x Tracheostomy Tube Holder
 - or
 - KACLIP® H** REF 3200x
1 x Tracheostomy Tube Holder
 - or
 - KACLIP® K** REF 3201x
1 x Tracheostomy Tube Holder
 - or
 - KACLIP® ClearH** REF 32015
1 x Tracheostomy Tube Holder

Product	Lenght	Width	REF
KACLIP® HM – beige	30 - 50 cm		32000
KACLIP® HM PED	13 - 25 cm		32001
KACLIP® H – beige	30 - 50 cm		32005
KACLIP® H – white	30 - 50 cm		32005-02
KACLIP® K – white	30 - 50 cm		32010
KACLIP® K – beige	30 - 50 cm		32011
KACLIP® Clear H	29 - 52 cm		32015

