

LARYNGOTEC® KOMBI STOMA BUTTON / FIX

- ▶ Stoma Button to keep the tracheostoma of laryngectomized patients open while minimizing the contact surface with the delicate tracheal mucosa
- ▶ For long-term use, can be used up to 3-6 months
- ▶ Medical grade silicone: soft, pliable, lightweight
- ▶ With 22 mm adapter (KOMBI) to connect to HMEs and/or speaking valves
- ▶ With (LARYNGOTEC® KOMBI Stoma Button FIX) or without (LARYNGOTEC® KOMBI Stoma Button) two lateral eyelets to fasten a tube holder to prevent dislocation of the stoma button
- ▶ Transparent, inconspicuous design
- ▶ Including broad adjustable OPTIFLAUSCH® tube holder (only variant „FIX“) and a product ID card



Order information:

LARYNGOTEC® KOMBI Stoma Button FIX, SIZE 8 - 12	REF 15950
LARYNGOTEC® KOMBI Stoma Button FIX short, SIZE 8 - 12	REF 15960
LARYNGOTEC® KOMBI Stoma Button, SIZE 10	REF 15955
LARYNGOTEC® KOMBI Stoma Button short, SIZE 11	REF 15965

When ordering, please state the size after the item number!

Size	Version	O.D. neck flange mm	I.D. mm	Length mm	REF FIX	REF
8	standard	11.0	8.0	22	15950-08	-
8	short	11.0	8.0	15	15960-08	-
9	standard	12.0	9.0	22	15950-09	-
9	short	12.0	9.0	15	15960-09	-
10	standard	13.0	10.0	22	15950-10	15955-10
10	short	13.0	10.0	15	15960-10	-
11	standard	14.0	11.0	22	15950-11	-
11	short	14.0	11.0	15	15960-11	15965-11
12	standard	15.0	12.0	22	15950-12	-
12	short	15.0	12.0	15	15950-12	-

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = length of shaft