


LARYVOX® EMERGENCY-SET TRACHEOSTOMY

Complete Set for the acute care of laryngectomized patients in case of emergency

- ▶ 1x LARYVOX® SOS safe-mask
Emergency ventilation mask for laryngectomees
- ▶ 1x TRACHEOTEC® VARIO CUFF
Tracheostomy tube with cuff, size 07
- ▶ 1x LARYVOX® Tracheospreizer (tracheal dilator)
Using the LARYVOX® Tracheospreizer, the cannula can be inserted safely and as quickly as possible
- ▶ 1x LARYVOX® Connect
Connects the 22 mm adapter (KOMBI) of the SOS safe-mask with a 15 mm connector (UNI) of e.g. a tracheostomy tube
- ▶ 5x LARYVOX® O₂ Connect
With multiadapter; to be connected to the Laryvox Tape base plate with 22 mm adapter (KOMBI) or to a tracheostomy tube with 15 mm (UNI) connector. The oxygen supply port allows the administration of extra oxygen
- ▶ 5x LARYVOX® Tape COMFORT / XL Oval
Base plate with strong adhesive strength and an extra-large surface, with 22 mm adapter (KOMBI)
- ▶ 3x OPTIFLUID®
Stoma oil towels with colourless neutral oil to lubricate the Tracheotec tracheostomy tube before insertion
- ▶ 1x Oxygen connection tube
- ▶ 1x Single-use syringe / 10ml



For all components of this kit separate and detailed data sheets are available.

Order information:

LARYVOX® EMERGENCY-SET Laryngectomy PRO

REF 48007