










# LARYVOX® HANDS-FREE VALVE LOW RESISTANCE / LARYVOX® HANDS-FREE VALVE KOMBI HME LOW RESISTANCE

- ▶ Transparent, lightweight
- ▶ Adjustable silicone membrane for hands-free speaking (laryngectomized patients with voice prosthesis)
- ▶ Very flexible thin silicone membrane facilitates low resistance
- ▶ Fits on 22 mm adapters (KOMBI)
- ▶ LARYVOX® Hands-Free Valve Low Resistance with handle to allow fine-tuning of the valve and to prevent clothes from occluding the valve
- ▶ LARYVOX® Hands-Free Valve Kombi HME Low Resistance with housing to allow the connection of 22 mm HMEs like e.g. LARYVOX® HMEs or HUMIDOTWIN®
- ▶ Individually packed



 Single patient use  
 Consult instructions for use  
 CE marking

## Order information:

**LARYVOX® Hands-Free Valve Low Resistance**

**LARYVOX® Hands-Free Valve Kombi HME Low Resistance**

**REF 21984**

**REF 21982**