

# LARYVOX® STARTER SET

## Contents:

- ▶ LARYVOX® TAPE  
5 x 3 pieces, different types, oval and rectangular shape, for different types of tracheostomas with different skin characteristics
- ▶ LARYVOX® HME  
15 pieces, low profile, inconspicuous, skin coloured, with stable ring housing for finger occlusion, with low respiratory resistance
- ▶ LARYVOX® HME HighFlow  
15 pieces, low profile, discreet design, skin coloured, with stable plastic housing for low pressure finger occlusion with low respiratory resistance
- ▶ LARYVOX® SECUTRACH® shower guard  
1 piece, with multiadaptor to connect to 15 mm connectors (UNI) (e.g. tracheostomy tubes) as well as to 22 mm adapters (KOMBI) (e.g. LARYVOX® TAPE)
- ▶ OPTIGARD® Skin Protection Film  
15 pieces, for preparative skin care when using Laryvox base plates, helps to protect the skin against irritations and increases the adhesive power of the base plate (e.g. LARYVOX® TAPE)
- ▶ OPTICLEAR® Adhesive Remover  
15 pieces, to quickly and gently remove residual skin adhesives after the removal of the base plate (e.g. LARYVOX® TAPE)



Please refer to the pictograms in the respective instructions of each product.

## Order information:

LARYVOX® Starter Set

REF 48810