



- ▶ To prevent ingress of water into the trachea during showering
- ▶ Opening for air supply tilted downward in a 90° angle to prevent the ingress of water while at the same time providing sufficient air to allow breathing without anxiety
- ▶ With multi-adapter: Can be connected to devices with 15 mm connectors (UNI), e.g. to tracheostomy tubes as well as to devices with 22 mm adapters (KOMBI), e.g. to base plates like LARYVOX® Tape Flexible
- ▶ Tight connection between tracheostomy tube or base plate and the LARYVOX® SECUTRACH®
- ▶ LARYVOX® SECUTRACH® must never be used while taking a full bath or during water sports activities



Single patient use only



Consult instructions for use



CE marking

Order information:

LARYVOX® SECUTRACH® shower guard

REF 47400