



OPTIFAHL® STOMA CLEANING WIPES

- ▶ Cleaning wipe for the cleaning and care of the sensitive skin around the tracheostoma
- ▶ Soaked in a mild cleaning solution, with natural and herbal ingredients
- ▶ Excessive skin oils, dirt particles, as well as remainders of skin adhesives are gently removed from the skin with this soft and lint-free wipe
- ▶ Also to be used for the preparative care of the skin when using protective LARYNGOFIX® cover foams (cf separate data sheet) or base plates such as LARYVOX® Tape (cf separate data sheet)
- ▶ Individually packaged in one-each sachets – ideal for travelling!
- ▶ Available in two different sizes 135 x 220 mm and 185 x 190 mm
- ▶ Also available in a dispenser for home-use:
The wipes can be removed easily and hygienically one-by-one from the dispenser; a lid prevents the wipes from drying out



Do not reuse



Consult instructions for use



CE marking

Order information:

OPTIFAHL® Stoma cleaning wipes, 185 x 190 mm, Pack of 30 (Individually packed)

REF 33200

OPTIFAHL® Stoma cleaning wipes, 135 x 220 mm, Pack of 30 (Individually packed)

REF 33210

OPTIFAHL® Stoma cleaning wipes, Container with 60 wipes

REF 33260