

OPTIFLAUSCH®

REF 32550



REF 32551



- ▶ Tracheostomy tube holder
 - "K" variant: with adherent fastener
 - "H" variant: with hook fastener
 - "K PED" variant: with adherent fastener, for pediatric use
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Broad and soft for high wearing comfort
- ▶ Approx. 4 mm thick and 2.5 cm wide.
- ▶ The tube holder is wrapped around the neck and the adherent fastener enables it to be easily and quickly adjusted to the desired neck size (2 parts, adjustable).
- ▶ Colour: white, beige
- ▶ Individually packed

SCOPE OF DELIVERY

- ▶ Manual

OPTIFLAUSCH® K **REF 3255x**
1 x Tracheostomy Tube Holder

or

OPTIFLAUSCH® H **REF 32551**
1 x Tracheostomy Tube Holder

Product	Length	Width	REF
OPTIFLAUSCH® K	24 - 45 cm	2,5 cm	32550
OPTIFLAUSCH® H	24 - 47 cm	2,5 cm	32551
OPTIFLAUSCH® K beige	24 - 45 cm	2,5 cm	32553
OPTIFLAUSCH® K PED	16 - 30 cm	2,5 cm	32555

