



# RESPITRACH® PERMANENT SUCTION CATHETER

- ▶ Closed Suction System for long-term tracheal suctioning
- ▶ For removing collected secretions from the lower respiratory passages
- ▶ To be used for ventilated patients – allows hygienic endotracheal aspiration of tracheal secretions without the need to disconnect the ventilation system from the tracheostomy tube/tracheal tube or to interrupt ventilation.
- ▶ Reduces the risk of infections (e.g. MRSA)
- ▶ Can be used up to 48 hours
- ▶ Catheter can be rinsed inside its closed sheath after with sterile saline solution each suctioning process
- ▶ Pack size of 8 each



## Order information:

**RESPITRACH® Permanent suction catheter SIZE 10-16**

When ordering, please state the size after the item number!

**REF 68400**