

(TRACHEO-)STOMA CATHETER METRIC S + L

25 + 50 cm



REF 68000/68900 CE 0482



REF 68100 CE 0482



REF 60700 CE 0482

- ▶ Sterile tracheostoma catheters for suctioning tracheal secretions
- ▶ Made of soft, flexible and transparent plastic
- ▶ With central opening and four lateral openings (“eyes”) at the tip to avoid adherence to the trachea
- ▶ With grading in 1 cm increments to facilitate control of the insertion depth
- ▶ Wide range of sizes (CH 6-18)
- ▶ Colour coded connectors
- ▶ Available in two different lengths: 25 cm and 50 cm (the short variant is particularly suited for tracheal suctioning)
- ▶ To be used with the TRACHFLOW® Fingertip, to control the suction power. The Fingertip is connected between the suction device and the sterile stoma catheter. One size of the TRACHFLOW® Fingertip fits all sizes of stoma catheters.

Article	Size	Length (cm)	Shore	Unit	Sterile	Atraumatic	REF
Tracheostoma catheter Metric S, 25 cm CH 8-16	8-16	25	A64	30	✓	✓	68100
Stoma catheter Metric L, 50 cm CH 6-18	6-18	50	A64	30	✓	✓	68000
Absaugkatheter Standard Metric L, 50 cm, CH 6-18	6-18	50	A78	100	✓	✓	68900

Order Information:

Tracheostoma catheter Metric S, 25 cm	CH 8-16	Pack of 30	REF 68100
Stoma catheter Metric L, 50 cm	CH 6-18	Pack of 30	REF 68000
Absaugkatheter Standard Metric L, 50 cm	CH 6-18	Pack of 100	REF 68900
TRACHFLOW® Fingertip		Pack of 10	REF 60700