

TRACHEOTEX® SCARF

- ▶ To cover and protect the tracheostoma of laryngectomized and tracheostomized patients
- ▶ To be tied around the neck using a broad, approx. 1 cm wide cotton strap
- ▶ A centrally incorporated insert, made of multiple layers of fine mesh warms, filters and humidifies the inspired air
- ▶ With adherent fastener
- ▶ Adjustable: An adjustment clip on the fixing strap allows the adjustment to different neck sizes (approx. 41 - 47 cm)
- ▶ Washable and reusable
- ▶ Available in many different colours and patterns matching the wardrobe (cf color sample card – available as extra brochure)



	Keep away from sunlight
	Keep dry
	Single patient use
	Consult instructions for use
	CE marking
	Hand wash
	Do not dry in the drum dryer

Order information:

TRACHEOTEX® Scarf, 41 - 47 cm

Please specify color number on the color sample card when ordering

REF 42000