


TRACHFLOW® METRIC S + L

25 + 50 cm



CE 0482



CE 0482

- ▶ Sterile, atraumatic tracheostoma catheters for particularly gentle suctioning of tracheal secretions
- ▶ With rounded tip to minimize the risk of tracheal injury
- ▶ Made of soft, flexible and transparent plastic
- ▶ With central opening and four lateral openings (“eyes”) at the tip to avoid adherence to the trachea
- ▶ With grading in 1 cm increments to facilitate control of the insertion depth
- ▶ Wide range of sizes (CH 8 -16)
- ▶ Colour coded connectors
- ▶ Available in two different lengths: 25 cm and 50 cm (the short variant is particularly suited for tracheal suctioning)
- ▶ To be used with the TRACHFLOW® Fingertip, to control the suction power.
The Fingertip is connected between the suction device and the sterile stoma catheter.
One size of the TRACHFLOW® Fingertip fits all sizes of stoma catheters.

Article	Size	Length (cm)	Shore	Unit	Sterile	Atraumatic	REF
TRACHFLOW® Metric S, 25 cm CH 10-16	8-16	25	A64	30	✓	✓	68630
TRACHFLOW® Metric L, 50 cm CH 10-16	8-16	50	A64	30	✓	✓	68730

Order Information:

TRACHFLOW® Metric S, 25 cm	CH 8 - 16	Pack of 30	REF 68630
TRACHFLOW® Metric L, 50 cm	CH 8 - 16	Pack of 30	REF 68730
TRACHFLOW® Fingertip		Pack of 10	REF 60700