

TRACHEOFIRST® COMPACT

SUCTION DEVICE

- ▶ Particularly lightweight suction device (approx. 2.2 kg) for domestic and hospital use for patients with normal levels of secretion production
- ▶ The electric swing-piston pump provides a suction power of 24 l/min and a maximum vacuum of up to 0.76 bar
- ▶ Including secretion container with a capacity of 1.000 ml
- ▶ Bacterial filter - prevents microorganisms and secretions from penetrating into the device
- ▶ Air filter - prevents microorganisms and secretions from being blown out of the device, and thus serves to protect the users
- ▶ The TRACHEOFIRST® COMPACT can also be operated via the 12 V car adapter in the car



Technical Data	
Suction power	24 l/min \pm 2 l/min
Max. vacuum	-76 kPa* (-760 mbar; -570 mmHg) \pm 4kPa
Secretion container	1,000 ml
Suction tube	Length: 1.80 m; diameter: 6 mm
Power supply	100-240 Vac (+/- 10%)

Technical Data	
Voltage frequency	50/60 Hz
Dimensions(HxWxD)	213 x 364 x 175 mm
Weight	approx. 2.2 kg

Order information:

TRACHEOFIRST® COMPACT

Bacterial filter for suction devices

TRACHFLOW® LINE PRO; Suction tube with fingertip, approx. length 1.30 m

Air filter pack of 10

REF 67600

REF 60800

REF 60506

REF 50000-07