




LARYVOX® EXTRA HME

- ▶ HME in discrete design with speech option via easy finger lock
- ▶ For use day and night (max. 24 hrs.)
- ▶ Excellent tactility: a rounded recess in the lid of the HME simplifies the locking mechanism and makes optimal positioning of the finger very easy to feel
- ▶ The lateral openings of the LARYVOX® EXTRA HME allow for unimpaired and easy breathing – even below clothing.
- ▶ A rabbet around the housing allows easy and quick removal of the HME (e. g. when coughing)
- ▶ With 22 mm adapter (KOMBI) to connect to any conventional base plates (e.g. LARYVOX® Tape) or silicone tubes for laryngectomees (e.g. LARYNGOTEC®)
- ▶ The following variants are available:

Normal - normal breathing resistance at normal physical exertion

Medium - medium breathing resistance at medium physical exertion

HighFlow - reduced breathing resistance at pronounced physical exertion



- ⊗ Do not reuse
- ☔ Keep dry
- ☀ Keep away from sunlight
- 📖 Consult instructions for use
- CE CE marking

Order information:

LARYVOX® EXTRA HME, Pack of 30

REF 49860

LARYVOX® EXTRA HME MEDIUM, Pack of 30

REF 49862

LARYVOX® EXTRA HME HighFlow, Pack of 30

REF 49861