



# LARYVOX® INSERTER EXPERT

Only for 16 Fr. / DUAL VALVE Shunt Valves


- ▶ The LARYVOX® INSERTER EXPERT consists of a placement instrument as well as an adapter
- ▶ LARYVOX® INSERTER EXPERT facilitates the perfect anterograde placement or replacement of a 16 Fr. or/and a DUAL VALVE shunt valve (16 Fr. and 20 Fr.)
- ▶ The placement instrument has three different fixation options to allow easy and safe placement of various types of shunt valves
- ▶ The applicator of the LARYVOX® INSERTER EXPERT enables effortless insertion of the shunt valve and is readily prepared in a few steps
- ▶ Additionally, the rounded applicator tip promotes gentle placement.



Applicator-Tip



Applicator-Opening for grasping of the placement instrument incl. shunt valve



Placement Instrument

Fixation Options



Consult instructions for use



Use by



Keep away from sunlight



Keep dry



CE marking



Do not reuse



Non-sterile

Product made in Germany

## Order information:

LARYVOX® INSERTER EXPERT

REF 25910