



# LARYVOX® INSERTER EXPERT SMALL

Only for 4-mm  
Shunt-Valves


- ▶ The LARYVOX® INSERTER EXPERT SMALL consists of a placement instrument as well as an adapter
- ▶ LARYVOX® INSERTER EXPERT SMALL facilitates the perfect anterograde placement or replacement of a 4 mm shunt valve
- ▶ The placement instrument has three different fixation options to allow easy and safe placement of various types of shunt valves
- ▶ The applicator of the LARYVOX® INSERTER EXPERT SMALL enables effortless insertion of the shunt valve and is readily prepared in a few steps
- ▶ Additionally, the rounded applicator tip promotes gentle placement




Consult instructions for use  
 Use by  
 Keep away from sunlight  
 Keep dry  
 CE marking  
 Do not reuse  
 Non-sterile  
**Product made in Germany**



Applicator-Tip



Applicator-Opening for grasping of the placement instrument incl. shunt valve



Placement Instrument

Fixation Options

## Order information:

**LARYVOX® INSERTER EXPERT SMALL**

**REF 25915**