


BLOM-SINGER® CLASSIC™ INDWELLING VOICE PROSTHESIS

- ▶ Blom-Singer® CLASSIC™ Indwelling is an indwelling voice prosthesis with integrated flap valve
- ▶ The prosthesis, which is made of medical-grade silicone, is furnished with a radiopaque ring at the oesophageal flange for a visual depiction
- ▶ Blom Singer® CLASSIC™ Indwelling Voice Prosthesis is available as 16 or 20 Fr. variant in size 4 - 20 (also available as sterile variant in 16 and 20 Fr., size 8 - 14)
- ▶ With the help of the LARYVOX® INSERTER the voice prosthesis can be placed anterograde quickly and easily:
 LARYVOX® INSERTER (REF 25900) - for different types of voice prosthesis
 LARYVOX® INSERTER EXPERT (REF 25910) - for 16 Fr. voice prosthesis (or DUAL VALVE™ 20 Fr.)
 LARYVOX® INSERTER EXPERT SMALL (REF 25915) - for 4 mm voice prosthesis

Application with gel capsule technique:

LARYVOX® INSERTER SPECIAL - (REF 25920) - for difficult tracheostoma



sterile



Single patient use



Consult instructions for use

Order Information:

| | | | |
|--|---------------|-----------|----------------------|
| Blom-Singer® CLASSIC™ Indwelling Voice Prosthesis | 16 Fr. | SIZE 4-20 | REF IN16XX-NS |
| | 20 Fr. | SIZE 4-20 | REF IN20XX-NS |
| Blom-Singer® CLASSIC™ Indwelling - sterile Voice Prosthesis | 16 Fr. | SIZE 8-14 | REF IN16XX-SO |
| | 20 Fr. | SIZE 8-14 | REF IN20XX-SO |
| LARYVOX® INSERTER | | | REF 25900 |
| LARYVOX® INSERTER EXPERT | | | REF 25910 |
| LARYVOX® INSERTER EXPERT SMALL | | | REF 25915 |
| LARYVOX® INSERTER SPECIAL | | | REF 25920 |