


BLOM-SINGER® PUNCTURE DILATOR TRACHEOESOPHAGEAL

- ▶ The Blom-Singer® Puncture Dilator Tracheoesophageal is a tracheoesophageal puncture dilator
- ▶ Fixed, conical silicone stent
- ▶ Used for successive expansion of the shunt canal or fistula canal while changing the prosthesis
- ▶ Prevents leakage and closure of the shunt channel or fistula channel when the shunt valve is removed
- ▶ The Blom-Singer® Puncture Dilator Tracheoesophageal is available with shaft thicknesses of 18 Fr. and 22 Fr. and is used for precise dilation for 16 and 20 Fr. voice prosthesis



Single patient use



Consult instructions for use

Order information:

Blom-Singer® Puncture Dilator Tracheoesophageal

18 Fr. (for 16 Fr. Voice Prosthesis)

BE6050

22 Fr. (for 20 Fr. Voice Prosthesis)

BE2050