

# TRACHEOPORT® JUNIOR

## SUCTION DEVICE FOR CHILDREN

- ▶ Suction device for domestic and hospital use for patients with normal levels of secretion production
- ▶ The optical appearance, which appeals to children, reduces dislikes / fears associated with inhalation
- ▶ The electric swing-piston pump provides a suction power of 24 l/min and a maximum vacuum of up to 0.76 bar
- ▶ Including secretion container with a capacity of 600 ml
- ▶ Bacterial filter - prevents microorganisms and secretions from penetrating into the device
- ▶ Air filter - prevents microorganisms and secretions from being blown out of the device, and thus serves to protect the users
- ▶ The TRACHEOPORT JUNIOR® suction device is supplied with a suitable storage bag (backpack) in addition to the complete set of accessories (consisting of power supply unit, secretion container, bacterial filter, air filters, suction tube, incl. clip, fingertip & car adapter)
- ▶ The device can be charged or operated with either the battery charger or the car charging cable
- ▶ The battery is fully charged in approx. 4 hours, operating time of battery approx. 60 minutes continuous operation
- ▶ When the battery charger or the car charging cable is connected, the full suction power is available, even if the battery is empty



Technical Data	
<b>Suction power</b>	24 l/min ± 2 l/min
<b>Max. vacuum</b>	-76 kPa* ( -760 mbar; -570 mmHg) ± 4kPa
<b>Secretion container</b>	600 ml
<b>Suction tube</b>	Length: 1.80 m; diameter: 6 mm
<b>Power supply (battery charger)</b>	100-240 Vac (+/- 10%)

Technical Data	
<b>Voltage frequency</b>	50/60 Hz
<b>Operating time with battery (without mains power supply)</b>	Operating time of battery approx. 60 min / continuous operation
<b>Charging time</b>	approx. 4 hrs.
<b>Dimensions (HxWxD)</b>	370 x 290 x 150 mm
<b>Weight</b>	approx. 2.4 kg



Observe instructions for use!



Device of protection class II



Typ BF applied part



Please observe the laws and regulations which apply to the disposal of electrical and electronic devices in the respective countries.



CE marking identification number of the notified body



IP classification

### Order Information:

**TRACHEOPORT® JUNIOR**

**Air filter, pack of 10**

**REF 63600**

**REF 50000-07**