








TRACHFLOW® SHUNT

10 cm

- ▶ Allows a targeted aspiration of shunt valves
- ▶ Made of soft, flexible and transparent plastic
- ▶ Large range (CH 8 and CH 10)
- ▶ Colour-coded connectors
- ▶ Length: 10 cm
- ▶ To be used in conjunction with the TRACHFLOW® Fingertip to control suction strength. The TRACHFLOW® Fingertip is used to connect the suction tube of the suction device to the sterile stoma catheter. The conical design allows one size to fit all sizes of stoma catheters.



 **0482** CE marking identification number of the notified body
 Do not reuse
 Do not use if package is damaged
 Keep dry
 Keep away from sunlight
 Sterilized using ethylene oxide

Order Information:

TRACHFLOW® SHUNT, CH 8, 10 cm

Pack of 30

REF 68800-08

TRACHFLOW® SHUNT, CH 10, 10 cm

Pack of 30

REF 68800-10

TRACHFLOW® Finger-Tip

Pack of 10

REF 60700