

# FAHL® NEBULIZER-SET WITH AEROSOLMASK

Inhalation



- ▶ The set consists of nebulizer and aerosol mask
- ▶ The set can be used at constant power, both vertically and at an angle
- ▶ The aerosol mask is easy to adjust
- ▶ FAHL® Nebulizer-Set with aerosolmask is a single-patient product and intended for single use
- ▶ Furthermore, the set is available in two different sizes (L and XL)
- ▶ Inhalation of steam, medications and other inhalants approved for use in inhalation therapy is possible
- ▶ The use for compressed air operated inhalation devices is possible

**Order information:**

**FAHL® NEBULIZER-SET WITH AEROSOLMASK**

Size L

**REF 51210-L**

**FAHL® NEBULIZER-SET WITH AEROSOLMASK**

Size XL

**REF 51220-XL**