


DURATWIX® JUNIOR NEO / NCF

- ▶ Tracheostomy tube for neonates, infants and toddlers (paediatrics), available without CUFF (REF 17432) and with CUFF (REF 17434)
- ▶ Adaptive, supportive and lightweight material
- ▶ Free of latex and DEHP
- ▶ Soft, flexible cannula tag with imprint of the essential specifications of the cannula
- ▶ Rounded tube tip to prevent irritation of the mucous membranes in the trachea
- ▶ With X-ray contrast strip
- ▶ Set includes broad a stable tube holder, an obturator and cannula passport



Order information:

DURATWIX® JUNIOR NEO

Size 2,5-4



REF 17432

DURATWIX® JUNIOR NCF

Size 2,5-4



REF 17434

When ordering, please state the size after the item number!

Size	O.D. mm	I.D. mm	Angel ca.	Cuff (Ø) mm	Length mm
2,5	4	2,5	120°	8	30
3	4,5	3	120°	9	31
3,5	5,1	3,5	120°	10	33
4	5,8	4	120°	12	35

O.D. = Outer diameter · I.D. = Inner diameter · Length = Length over middle curve