


DURATWIX® JUNIOR PED / PDC

- ▶ Tracheostomy tube for children (paediatrics), available without CUFF (REF 17436) and with CUFF (REF 17438)
- ▶ Adaptive, supportive and lightweight material
- ▶ Free of latex and DEHP
- ▶ Soft, flexible cannula tag with imprint of the essential specifications of the cannula
- ▶ Rounded tube tip to prevent irritation of the mucous membranes in the trachea
- ▶ With X-ray contrast strip
- ▶ Set includes broad a stable tube holder, an obturator and cannula passport



Order information:

DURATWIX® JUNIOR PED Size 2,5-6  **REF 17436**

DURATWIX® JUNIOR PDC Size 2,5-6  **REF 17438**

When ordering, please state the size after the item number!

Size	O.D. mm	I.D. mm	Angel ca.	Cuff (Ø) mm	Length mm
2,5	4,0	2,5	120°	8	32
3	4,5	3,0	120°	9	35
3,5	5,1	3,5	120°	10	39
4	5,8	4,0	120°	12	43
4,5	6,5	4,5	120°	12	47
5	7,1	5,0	120°	15	49
5,5	7,6	5,5	120°	15	54
6	8,3	6,0	120°	18	60

O.D. = Outer diameter · I.D. = Inner diameter · Length = Length over middle curve