

## FAHL® EMERGENCY-SET PED

Complete set for the emergency care of children:

- ▶ 1x LARYVOX® SOS SAFE-MASK  
Emergency ventilation mask for laryngectomees
- ▶ 1x DURATWIX® JUNIOR PED  
Tracheostomy tube for children, size 05
- ▶ 1x LARYVOX® TRACHEOSPREIZER PED  
Using the Tracheospreizer, the cannula can be inserted safely and as quickly as possible
- ▶ 5x LARYVOX® O2 HME  
HME filter with oxygen connection and multi-adapter for connecting 15 mm (UNI) connectors and 22 mm (KOMBI) Adapters
- ▶ 2x SENSOTRACH® DUO PED SLIT  
Tracheal compress serves to absorb tracheal secretions and to cushion the tracheal cannula shield against the skin
- ▶ 3x OPTIFLUID® STOMA OIL TOWEL  
Stoma oil towels with colourless neutral oil to lubricate the Tracheotec tracheostomy tube before insertion
- ▶ 3x OPTIFAHL® STOMA CLEANING WIPES  
Stoma cleaning wipes for cleaning and caring for the skin area around the tracheostoma
- ▶ 1x Oxygen connection tube
- ▶ 1x Single-use syringe / 10ml
- ▶ 1x FAHL® CATHETER MOUNT  
Catheter mount (length 15 cm) for the connection between tracheal cannula and ventilation system
- ▶ 3x Tracheostoma Catheter Metric S  
Tracheostoma suction catheter for gentle aspiration of tracheal secretions



For all components of this kit separate and detailed data sheets are available.



### Order information:

**FAHL® EMERGENCY-SET PED**

**REF 48009**